# HELMUT ZEPF

[www.zepf-dental.com](http://www.zepf-dental.com)

Close to you.

MANUFACTURER SINCE 1921

[www.zepf-dental.com/  
anwenderinformationen](http://www.zepf-dental.com/anwenderinformationen)

CE 0483

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## Reprocessing of resterilizable ZEPF hand-held instruments and their accessories

### GENERAL REMARKS

Instruments may only be used by persons who are specially trained for this purpose. This also applies to reconditioning instruments. According to the brochure of the Working Group Instrument Reprocessing (AKI, yellow and red brochures), new instruments always have to pass through the complete cycle of instrument reprocessing at least once prior to first use, as described in the section “**Steps in Instrument Reprocessing**”. All instruments must be cleaned, disinfected and sterilized before each use. Following removal from the protective packaging, all non-sterile instruments must also be cleaned, disinfected and sterilized before first use. Effective cleaning and disinfection constitute an absolute precondition for proper sterilization of the instruments.

The user is responsible for the sterility of the instruments. Therefore, make sure only validated methods are used for cleaning, disinfection and sterilization. Sterilization equipment must also be serviced and checked at regular intervals. The validated parameters for the cleaning and sterilization cycles must also be checked regularly.

Also comply with the legal requirements valid in your country and the hygiene regulations of the medical practice or hospital. Improper handling and care or use for purposes for which they are not intended may result in premature wearing out of surgical/dental instruments. Persons using these instruments should be familiar with use and handling of surgical/dental instruments, accessories and related devices.

### INSPECTION AND FUNCTIONAL CHECK

It is very important to inspect each surgical/dental instrument for ruptures, cracks or malfunction prior to each use. Above all, areas such as the cutting edges, tips, closures, locks and catches, as well as all moving parts, must be carefully checked. Never use damaged instruments. Do not carry out repairs yourself. Service and repairs must be done by properly trained personnel only. If you have any questions in this regard, ask the manufacturer or your technical medical department.

### INSTRUMENT TRAYS

The **HELMUT ZEPF Protector System** is a proven approach with considerable advantages. It provides an ideal method for ready instrument layout according to the process involved (diagnostics, prophylaxis, extraction, etc.). The Tray-in-Tray System is an optimized approach to efficient cleaning, disinfection, sterilization and storage.

### STEPS IN INSTRUMENT REPROCESSING

#### RKI CLASSIFICATION

**Semi-critical:** Instruments that come into contact with mucosa or pathologically changed skin.

**Detail classification semi-critical “A”:** This group includes e.g. orthodontic forceps and other orthodontic instruments and/or instruments intended solely for this application.

**Detail classification semi-critical “B”:** Rotating instruments not intended for invasive surgical applications. This group includes, e.g.: Instruments for general preventive, restorative or orthodontic treatment.

**Critical medical devices:** Invasive surgical devices. Medical devices for use of blood, blood products and other sterile medicinal products and medical devices that penetrate skin or mucosa and thereby come into contact with blood, internal tissues or organs, including wounds.

**Detail classification critical “A”:** No special reprocessing specifications. This group includes e.g. universal procedures, probes, PA probes, tweezers, PA curettes, chisels, hoes, PA scalpels, endo instruments, forceps, clamps and needle holders, filling and modeling instruments, excavators, matrix retainers, Benex extractors, root extractors, periostomes, spacers and retractors, scissors, raspatories, scrapers, osteotomes, bone mills, sinus lift elevators and raspatories, measuring instruments.

**Detail classification critical “B”:** With more stringent reprocessing specifications since cleaning effectiveness cannot be directly assessed by inspection. This group includes e.g. medical suction devices, water syringes, hollow-cylinder osteotomes, trephines and membrane punches.

### PREPARATION AND TRANSPORT

The products must be stored and transported to the reprocessing site in a closed container to avoid damage to the instruments and contamination of the environment.

### CLEANING AND DISINFECTION

#### Principles

For cleaning and disinfection, if possible an automated procedure (washer-disinfector or WD) should be used. A manual procedure – even using an ultrasonic bath – should only be used if an automated procedure is not available or in accordance with country-specific requirements (e.g. in Germany the automated procedure is mandatory for critical B-products) due to the significantly lower effectiveness and reproducibility.

#### Pre-treatment

Coarse impurities must be removed from the products immediately after use (within 2 h maximum). If, due to the duration of use or as a result of organizational aspects, this time cannot be adhered to, the user must determine and validate measures autonomously to prevent the contamination from drying on:

1. Where applicable, dismantle the products as far as possible.
2. Rinse the products at least 1 min. under running water (temperature < 35°C/95°F). Moving parts should be shaken back and forth at least five times during pre-rinsing. Where applicable, brush and rinse all lumens of the products at least five times. Open and close the joints at least five times.
3. Place the disassembled products in a sufficiently large pre-cleaning bath (in an ultrasonic bath that has not been activated yet) for the specified soaking time so the products are completely submerged. Ensure that the products do not touch each other. Support pre-cleaning by completely brushing all internal and external surfaces at the beginning of the soaking time. The joints are in the open position.

#### Caution

Be careful with products with narrow gaps in which brush bristles can get stuck! Brushes for the channels must be slightly larger than the respective channel inside diameter; the shaft length of the brush must be at least as long as the channel. Moving parts should be shaken back and forth at least five times during pre-cleaning. If applicable: Rinse all lumens of the products at least five times at the beginning and end of the soaking time (aids and minimum volume depending on the cavity to be rinsed).

4. Activate the ultrasound for another minimum soaking time (but not less than 5 min.).
5. Then remove the products from the pre-cleaning bath and rinse thoroughly (at least 1 min.) with water at least five times. Moving parts should be shaken back and forth at least five times during post-rinsing. If applicable: Rinse all lumens of the products at least five times (aids and minimum volume depend on the cavity to be rinsed).

When selecting the cleaning agent<sup>2</sup> to be used, please take into consideration that

- it is generally suitable for cleaning invasive medical devices made of metals and plastics
- the cleaning agent is suitable for ultrasonic cleaning (no foam formation)
- the cleaning agent is compatible with the products (see chapter “Material Stability”).

The concentrations, temperatures and soaking times, as well as rinsing instructions specified by the manufacturer of the cleaning agent or cleaning and disinfecting agent, must be strictly adhered to. Only use freshly prepared solutions, only sterile or low-bacteria

## Reprocessing of resterilizable ZEPF hand-held instruments and their accessories

(max. 10 bacteria/ml) and low endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water<sup>1</sup>) and only a soft, clean and lint-free cloth and/or filtered air.

### Caution

Take care in case of products with rough surfaces, threads, sharp edges or comparable aspects on which particles from the cloth can be caught!

### AUTOMATED CLEANING / DISINFECTION (washer-disinfector or WD)

When selecting the WD, ensure that

- the WD generally has verified effectiveness (e.g. the German Society for Hygiene and Microbiology (DGHM) or FDA approval/clearance/registration or CE marking according to DIN EN ISO 15883)
- if possible, a tested program for thermal disinfection (A0 value > 3000 or – for older devices – at least 5 min. at 90°C/194°F) is used (in chemical disinfection, danger of disinfecting agent residues on the products)
- the program used is suitable for the products and includes sufficient rinsing steps (at least three degrading steps after cleaning (or neutralization if applied) or conductance-based rinsing control recommended to effectively prevent detergent residues)
- only sterile or low-bacteria (max. 10 bacteria/ml) and low endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water)<sup>1</sup> is used
- the air used for drying is filtered (oil-free, low-bacteria and low-particle) and
- the WD is regularly maintained, checked and calibrated.

When selecting the cleaning system to be used, ensure that

- it is generally suitable for cleaning invasive medical devices made of metals and plastics
- providing no thermal disinfection is used – a suitable disinfecting agent with verified effectiveness (e.g. the German Association for Applied Hygiene (VAH) / The German Society for Hygiene and Microbiology (DGHM) or FDA/EPA approval/clearance/registration or CE label) is also used and that this is compatible with the cleaning agent used and
- that the chemicals used are compatible with the products (see chapter “Material Stability”).

The concentrations, temperatures and soaking times, as well as rinsing instructions specified by the manufacturer of the disinfecting agent, must be strictly adhered to.

### Procedure:

1. Place the disassembled products coming from the pre-cleaning into the WD. Ensure that the products do not touch each other, joints are open and, if necessary, are connected to the rinsing connection.
2. Start the program.
3. Remove the products from the WD at the end of the program.
4. Inspect and pack the products as soon as possible after removal (see chapter “Inspection” and “Packaging”, possibly after additional drying in a clean area).

Verification of the products’ general suitability for effective automated cleaning and disinfecting was provided by an independent, officially accredited and recognized (§ 15 (5) German Medical Devices Act (MPG)) test laboratory using the WD G 7836 CD (thermal disinfection, Miele & Cie. KG, Gütersloh) and the neodisher MediClean forte pre-cleaning and cleaning agent (Dr. Weigert GmbH & Co. KG, Hamburg). Here, the procedure described above was taken into consideration.

### MANUAL CLEANING AND DISINFECTION

When selecting the cleaning and disinfecting agent to be used, ensure that

- it is generally suitable for cleaning invasive medical devices made of metals and plastics

- the cleaning agent is suitable for ultrasonic cleaning (no foam formation)
- a suitable disinfecting agent with verified effectiveness (e.g. The German Association for Applied Hygiene (VAH) / The German Society for Hygiene and Microbiology (DGHM) or FDA/EPA approval/clearance/registration or CE marking) is also used and that this is compatible with the cleaning agent used and
- that the chemicals used are compatible with the products (see chapter “Material Stability”).

Combined cleaning/disinfecting agents should not be used if possible.

The concentrations, temperatures and soaking times specified by the manufacturer of the cleaning and disinfecting agent, must be strictly adhered to. Only use freshly prepared solutions, only sterile or low-bacteria (max. 10 bacteria/ml) and low endotoxin (max. 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water)<sup>1</sup> and only a soft, clean and lint-free cloth and/or filtered air.

### Caution

Be careful with products with rough surfaces, sharp edges or comparable aspects on which particles from the cloth can be caught!

### Procedure: Cleaning

1. Place the disassembled products coming from pre-cleaning in a sufficiently large cleaning bath (in an ultrasonic bath that has not been activated yet) for the specified soaking time so that the products are completely submerged. Ensure that the products do not touch each other. Support the cleaning process by completely brushing all internal and external surfaces with a soft plastic brush and shake moving parts back and forth several times. The joints are in the open position.
2. Activate the ultrasound for another minimum soaking time (but not less than 5 min.).
3. Then remove the products from the cleaning bath and rinse thoroughly (for at least 1 min.) with water at least five times. Moving parts should be shaken back and forth several times during post-rinsing.
4. Inspect the products (see chapter “Inspection”).

### Disinfection

5. Place the products for the predefined soaking time in the cleaning bath so that the products are completely submerged. Ensure that the products do not touch each other. The joints are in the open position.
6. Then remove the products from the disinfection bath and rinse thoroughly (for at least 1 min.) with water at least five times. Move moving parts back and forth several times during rinsing.
7. Dry the products by blowing off/out with filtered compressed air.
8. Pack the products as soon as possible after removal (see chapter “Packaging”, possibly after additional drying in a clean area).

### Caution

During manual cleaning and disinfection, where there is a potential risk of injury and infection, additional occupational safety measures (e.g., protective clothing, safety goggles, gloves, air filtration) must be observed in accordance with national regulations (e.g., TRBA [Technical Rules for Biological Agents] 250 in Germany).

### Inspection

Check all products after cleaning or cleaning/disinfection for corrosion, damaged surfaces, chippings, contaminants and stains and discard damaged products. Products that are still contaminated must be cleaned again and disinfected.

### Maintenance

Reassemble disassembled products.

In the case of oiling joints, care must be taken to ensure that only instrument oils (white oil, without further additives) are used which – taking into consideration the maximum sterilization temperature applied – are approved for steam sterilization and have tested biocompatibility, and that only a small amount is applied to the joints.

## Reprocessing of resterilizable ZEPF hand-held instruments and their accessories

CE 0483

### FUNCTIONAL TESTING

Check articulated instruments for ease of movement (avoid excessive play). The function of locking mechanisms (locking wheel) should be checked. **All instruments:** Perform visual inspection for damage and wear. Cutting edges should not have nicks and should be even. Check long, slender instruments (especially articulated instruments) for damage. If instruments are part of a larger construction, this must be checked together with the associated components.

### Packaging

If necessary, sort the cleaned and disinfected products into the corresponding sterilization tray. Please package the products or the sterilization trays in sterilization containers, or very large products in single-use sterilization packaging (single or double packaging) that meet the following requirements (material/process):

- DIN EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature stability up to at least 138°C (280°F) sufficient steam permeability)
- sufficient to protect the products or sterilization packaging from mechanical damage
- regular maintenance according to the manufacturer's specifications (sterilization container)
- do not exceed a maximum weight of 10 kg per pack / content of the sterilization container.

### STERILIZATION

For sterilization, only the following sterilization methods may be used; other sterilization methods are not allowed.

#### Steam sterilization

- Fractionated vacuum procedure<sup>3,4</sup> (with sufficient product drying<sup>5</sup>)
- Steam sterilizer in accordance with DIN EN 13060/DIN EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- Validation in accordance with DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance assessment (PQ))
- Maximum sterilization temperature 134°C (273°F); plus tolerance in accordance with DIN EN ISO/ANSI AAMI ISO 17665
- Sterilization time (soaking time at sterilization temperature):

| Country         | fractional vacuum procedure                                                                                    |
|-----------------|----------------------------------------------------------------------------------------------------------------|
| Germany         | at least 5 min. <sup>6</sup> at 134°C (273°F)                                                                  |
| USA             | min. 4 min. at 132°C (270°F), drying time min. 20 min <sup>4</sup>                                             |
| France          | at least 5 min. <sup>6</sup> at 134°C (273°F)<br>if required for prion inactivation sterilization time 18 min. |
| other countries | at least 5 min. <sup>6</sup> at 132°C (270°F) / 134°C (273°F)                                                  |

Gravity displacement is not recommended<sup>4</sup>.

Verification of the general suitability of the products for effective steam sterilization has been provided by an independent, officially accredited and recognized (§ 15 (5) MPG) test laboratory using the HST 6x6x6 steam sterilizer (Zirbus technology GmbH, Bad Grund) and using the fractionated vacuum procedure and LAWTON MEDOIL (oils for joints and friction surfaces). Here, typical conditions in clinics and medical practices and the procedure described above were taken into consideration.

The flash sterilization process is generally not permitted.

Do not use hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, or plasma sterilization.

### Material Stability

When selecting cleaning and disinfecting agents, please ensure that they do not contain the following components:

- organic, mineral and oxidizing acids (minimum permissible pH value 5.5)
- Alkalis/strong alkaline solutions (neutral/enzymatic max. permissible pH 8.5, mandatory requirement for products made of aluminum, Ferrozell or other alkali-sensitive materials) or alkaline cleaner (max. permitted pH 11, mandatory requirement for products with intended use in prion-critical areas, e.g. in accordance with Appendix 7 of the KRINKO RKI BfArM recommendation for reprocessing)
- Organic solvents (e.g. alcohols, ethers, ketones, benzines)
- Oxidizing agents (e.g. hydrogen peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic/halogenated hydrocarbons

Never clean products, sterilization trays or sterilization containers with metal brushes or steel wool.

All products, sterilization trays and sterilization containers may only be exposed to temperatures under 138°C (280°F).

### STORAGE

After sterilization, the products must be stored dry and free of dust in the sterilization packaging at room temperature. Avoid large temperature fluctuations. This prevents the formation of condensate and the related corrosion.

In the case of proper storage, there is no shelf-life specification or functional restriction for instruments supplied non-sterile after production.

### EXCEPTIONS!

#### See the reprocessing specifications for the items in the

**following list:** The instruments listed below are made in part from chrome-plated parts for technical reasons and must not be placed in the thermal disinfectant or ultrasound bath.

- Cylinder ampoule syringe thumb ring
- Replaceable ampoule holder
- Ointment cannulas
- Napkin chain
- Mouth mirror
- Resection mirror
- Chrome-plated instruments

### CUSTOM-MADE INSTRUMENTS

HELMUT ZEPF Medizintechnik GmbH develops and manufactures instruments according to customer specifications to produce custom-made products that satisfy individual customer needs. If not specific information is provided, maintenance and care of these products must also follow these instructions.

### GUARANTEE

HELMUT ZEPF Medizintechnik GmbH supplies only tested and defect-free products to its customers. All of our products are designed and manufactured to comply with the most stringent quality requirements. We assume no liability for products that have been altered compared to the original, used for purposes for which they were not intended or used improperly. Products that are repaired or altered at any company other than HELMUT ZEPF or a company commissioned for this purpose by HELMUT ZEPF are no longer authorized to bear the CE symbol.

## Reprocessing of resterilizable ZEPF hand-held instruments and their accessories

CE 0483

**Note regarding return shipments:** Due to legal regulations and, above all, to protect our employees, a signed proof of decontamination is required for each return shipment. This form must be completed, signed and attached to all product returns (complaints / repairs / other reasons for return). Please pack the products in a way to avoid any risk of injury for our staff when unpacking.

The form is available for download on our website  
„[https://www.zepf-dental.com/IFU/080300008\\_GB.pdf](https://www.zepf-dental.com/IFU/080300008_GB.pdf)“.

### Caution

All serious incidents occurring related to the product must be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is resident.

The above instructions have been validated as SUITABLE for preparation of a medical device for its reuse by the medical device manufacturer.

The reprocessor bears responsibility for achieving the desired results of the actual reprocessing with the equipment, materials and personnel used in the reprocessing facility for that purpose. This normally requires validation and routine monitoring of the process. Also, any deviation from the instructions provided on the part of the reprocessor should be carefully assessed for effectiveness and possible negative effects.

### FURTHER RECOMMENDATIONS / REFERENCES

- Ger. Federal Health Gazette (Bundesgesundheitsblatt)  
“Hygienic Specification for Reprocessing of Medical Devices”
- The present form can be downloaded from our website  
[www.zepf-dental.com](http://www.zepf-dental.com).

<sup>1</sup> Should you consider a lower water quality as sufficient against the background of national recommendations (e.g. in Germany KRINKO/RKI/BfArM recommendation for treatment), this is your sole responsibility.

<sup>2</sup> If you use a cleaning and disinfecting agent for this purpose - e.g. for occupational safety reasons - please note that it should be aldehyde-free (otherwise fixation of blood contamination), have a tested effectiveness (e.g. VAH/DGHM or FDA/EPA approval/cleanance/registration or CE marking), be suitable for disinfecting the products and be compatible with the products (see chapter “Material Stability”). Please note that the disinfecting agent used during pre-treatment is for personal protection only and cannot replace the later disinfection step which must be carried out after cleaning.

<sup>3</sup> at least three vacuum steps

<sup>4</sup> The use of the less effective gravity displacement is only permissible if the fractionated vacuum procedure is not available, it requires significantly longer sterilization times and must be validated by the user for specific products, devices, processes and parameters.

<sup>5</sup> The actual required drying time depends directly on parameters that are the sole responsibility of the user (loading configuration and density, sterilizer status, ...) and must, therefore, be determined by the user. Nevertheless, drying times should not be under 20 min.

<sup>6</sup> or extended sterilization time (e.g. 18 min.) for prion inactivation according to national regulations (not relevant for USA)

### Explanation of symbols

|                                                                                     |                                             |                                                                                       |                          |
|-------------------------------------------------------------------------------------|---------------------------------------------|---------------------------------------------------------------------------------------|--------------------------|
|  | Observe the electronic instructions for use |  | Unique device identifier |
|  | Manufacturer                                |  | Quantity                 |
|  | Date of manufacture                         |  | HIBC Barcode             |
|  | Catalog number                              |  | CE marking               |
|  | LOT no. (batch no.)                         |  | Caution                  |
|  | Non-sterile                                 |  | Medical device           |

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